



Overview and Technical Tips



www.ptglab.com



CONTENTS



- **3–7** Introduction
- 8–14 Transfection Techniques Cationic-lipid mediated transfection Calcium phosphate transfection Viral transfection Electroporation
- **15–16** FAQs and Troubleshooting
 - 17 Contact Us







- A transfection is a process that forces nucleic acids into a cell.
- Different methods are used for a transfection (chemical, physical or biological method).
- The transfected nucleic acids may stay transient or stable in the cell, dependent on the method of transfection.

Transfection:

- Nucleic acid transfer into eukaryotic cells.

Transformation:

- Non-viral DNA transfer into non eukaryotc cells, bacteria, plants.

Transduction:

- DNA transfer mediated by a virus.



Gene Expression Vs. Gene Inhibition

Gene Expression

- Transfection is used to overexpress a protein-of-interest.

Gene Inhibition

- Transfection is used to down-regulate a protein-of-interest via RNA interference (RNAi).
- miRNA, shRNA or siRNA can induce gene silencing.





Transient Vs. Stable Transfection

Transient Transfection

- The transfected nucleic acid stays in the cell for a limited period of time.
- The transfected material gets lost during cell division.
- Transfection effects can be observed normally between 1-7 days after transfection.
- Quick and easy.

Stable Transfection

- The transfected material is integrated into the genome.
- Uptake of foreigner DNA material into the genome requires optimal transfection conditions and selection of positive cells.
- Experiment needs approx. 1 month.







Main Factors Influencing A Transfection

Factor	Suggestion	
Cell health	Cell and transfection technique appropriate media should be used.	
Cell culture	Correct cell density is needed (40-70%).	
	Low cell passage.	
	Transfection efficiency might change with number of cell passage.	
Transfection material	Optimal amount of transfection material has to be validated for each experiment (cell type, cell number, transfection method).	
	Use high quality/pure material.	





In Conclusion

- These factors have to be considered before starting a transfection.



Overview of the most common transfection techniques			
Technique	+	-	
Cationic-lipid mediated transfection	Fast. Easy. Commercially available. Transient/stable transfection. Also for <i>in vivo</i> use suitable.	Some cell lines can't be transfected. Mainly serum sensitive. Optimization necessary, transfection efficiency highly depends on the cell line.	
Calcium phosphate transfection	Fast. Easy. Transient/stable transfection. Also for <i>in vivo</i> use suitable.	Transfection efficiency depends on cell type. Careful preparation of reagents. Cytotoxic. Reproducibility. Not for <i>in vivo</i> use.	



Overview of the most common transfection techniques		
Technique	+	-
Viral transfection	For difficult to transfect cells. High transfection efficiency.	Experienced user. Cytotoxic for cells. Biosaftly 2 lab might be required. Laborious experiment. Expensive.
Electroporation	Suitable for primary cells. Quite easy handling. Reproducible results. Fast set-up.	Initial acquisition costs. Optimization for each condition required.

10



Cationic-lipid Mediated Transfection The method relies on the formation of a complex between the DNA and a cationic lipid-based reagent. The complex is taken up by the cell via endocytosis.

- Fast and easy.
- Prize range of transfection reagents varies.
- Not working in all cell types.
- Transfection results may vary.
- Transfection success depends on quality of transfection reagent.





Calcium Phosphate Transfection

The methods is based on the formation of precipitates of calcium phosphate and DNA. The precipites are taken up by the cell via endocytosis.

- Easy protocol.
- Not expensive.
- Low transfection efficiency compared to other transfection reagents.
- Works well with most cell lines.
- Toxic for most primary cells.
- For transient and stable transfection.
- Sensitive to small changes in pH/ temperature, does not work with high phosphate buffers.





Viral Mediated Transfection (Transduction)

- Adeno-, Retro- or Lentivirus are used. _
- Alternative for difficult to transfect cells. _
- Suitable for protein overexpression _ or protein silencing.
- Additional pre-transduction steps _ are required:

Design of the virus mediating the transgene.

Production/purification of virus particles.

Virus titer estimation.





www.ptglab.com



Electroporation

The method uses an electric field that transiently permeabilizates the membrane. During this period the foreigner material can enter the cell.

- Works well with almost all cell types.
- High transfection efficiency.
- Requires optimization as cell permeability can lead to cell damage.
- Experienced user needed.
- Quite expensive method.







FAQS AND TROUBLESHOOTING



In Conclusion

- These factors have to be considered before starting a transfection.



FAQS AND TROUBLESHOOTING

Main Factors Influencing A Transfection

Factor	Suggestion	
Cell health	Cell and transfection technique appropriate media should be used.	
Cell culture	Correct cell density is needed (40-70%).	
	Low cell passage.	
	Transfection efficiency might change with number of cell passage.	
Transfection material	Optimal amount of transfection material has to be validated for each experiment (cell type, cell number, transfection method).	
	Use high quality/pure material.	



CONTACT US

Proteintech Group US Head Office proteintech@ptglab.com

Proteintech Europe United Kingdom europe@ptglab.com

Proteintech *China Office* service@ptglab.com

Support

Available 24 hours via Live Chat and 9–5 (CDT) via phone.

Please visit us at www.ptglab.com for more information about our antibodies and technical tips.



www.ptglab.com