

CHROMATIN IMMUNOPRECIPITATION

www.ptglab.com

Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

Day 1

ChIP & Cell Line Samples:

For suspension cells, ensure the appropriate number are resuspended in 10 ml of fresh media to which 37% of formaldehyde solution is subsequently added until a final concentration of 1% is reached. Vortex for a few seconds to displace from the bottom and incubate for 10 minutes at room temperature in agitation.

ChIP & Tissue Samples:

Grind frozen tissue into powder with a pestle and mortar. Pour the resulting powder into a 15 or 50 ml falcon. Add 10 ml PBS + 270 µl formaldehyde 37% (final concentration 1%) to the frozen powder. Vortex for a few seconds to displace from the bottom, next incubate for 10 minutes at room temperature in agitation.

- Block the reaction with 500 µl Glycine 2.5 M (final concentration 0.125 M). Incubate for 5 minutes at room temperature.
- Transfer the cells to a 50ml falcon and centrifuge at 2500 rpm for 5 minutes at 4°C.
- Discard the supernatant and wash twice with ice-cold **PBS** ph 7.4 and centrifuge at 2500 rpm for 5 minutes at 4°C after each washing.
- Resuspend cells in 5 ml **Cell Lysis Buffer** supplemented with protease inhibitor. To facilitate the cell membrane breaking, pass the lysate 3 times to a douncer. Incubate for 15 minutes at 4°C.
- Centrifuge at 4000 rpm for 5 minutes at 4°C. Discard the supernatant.
- Resuspend nuclei in **Cell Lysis Buffer**.

Tip 1

Pipette with cut tips to homogenize better.

- Divide the sample into small aliquots and sonicate for 15 minutes (high power; 30 seconds sonication, 30 seconds rest).

Tip 2

Put ice into the sonicator to avoid sample overheat.

- Centrifuge at 12000 rpm for 10 minutes at room temperature to remove nuclear debris. Discard the pellet. Repeat this passage the pellet cannot be detected. Store the samples at -20°C.

Tip 3

If SDS precipitates, dissolve it prior to centrifuge.

Tip 4

Take an aliquot of chromatin to quantify (2–3 µl) and to assess the size (30–40 µl).

CHROMATIN IMMUNOPRECIPITATION

www.ptglab.com

Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

Day 1

DNA Fragment Testing:

- De-crosslink chromatin by incubating samples at 65°C for 4 hours (results may improve with an overnight incubation).
- Incubate for 30 minutes with Proteinase K 50 µg/ml final concentration at 42°C.
- Add 1 volume of Phenol:Chloroform:Isoamyl Alcohol (25:24:1). Mix with vortex and let samples stand at room temperature for 2–3 minutes.
- Centrifuge at 12000 rpm for 5 minutes and transfer the aqueous phase to a new tube.

Tip 5

If the interphase is dirty, repeat steps c and d.

- Add 1/5 volume of AcNH₄ 10 M and 2.5 volumes of EtOH 100%. Mix and let the DNA precipitate for at least 30 minutes at -20°C.
- Centrifuge at max speed for 15 minutes at 4°C. Discard the supernatant.
- Wash with 70% ethanol and centrifuge at max speed for 15 minutes at 4°C. (Try to discard as much supernatant as you can without touching the pellet.)
- Resuspend with **TE Buffer** and pipette until complete dissolution.
- Incubate for 30 minutes at 37°C with RNase A at a final concentration of 50 µg/µl.
- Prepare a 1.5% agarose gel.

Day 2

- Take 70 µl of Magnetic Beads for each sample to be immunoprecipitated. (Take the extra volume in excess: 0.5–1 times more.)
- After precipitation with a magnet, discard the supernatant and wash twice with 600 µl at 5% BSA/PBS.
- After the second wash reconstitute the initial volume (70*N° of samples µl). (Take the extra volume in excess: 0.5–1 times more.)
- Take 20 µg (dependant on the tissue/cell type) for each sample, dilute the chromatin 1:10, and bring to a final volume of 1 ml with **Dilution Buffer**.
- Take 25 µl of beads for each sample and add them to the chromatin for the pre-clearing step.
- Divide the remaining beads into 45 µl aliquots. Add the corresponding antibody to each tube, plus a negative control (specific IgG). Incubate overnight at 4°C in a rotating wheel.

CHROMATIN IMMUNOPRECIPITATION

www.ptglab.com

Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

Day 3

Please Note: From this stage it is better to work with siliconized tubes.

- a. Discard the beads from the chromatin samples by putting the tubes in the magnet.
- b. Wash the Ab-Bead complexes twice with ice-cold 300 µl BSA/PBS 5%. Spin after second wash to displace the sample from the tube walls. Remove as much of the supernatant as possible.
- c. Add 1 ml of the chromatin to each sample and resuspend with the tip. Incubate for 2 hours at 4°C in a rotating wheel.

Tip 6

Store the excess chromatin from the input sample.

- d. Spin the samples and put them on the magnet.
- e. Wash twice with 1 ml **Low Salt Buffer**.
- f. Wash twice with 1 ml **High Salt Buffer**.
- g. Wash twice with 1 ml **LiCl Buffer**.
- h. Wash twice with 1 ml **TE** when adding the second washes. Change tubes for new ones.
- i. Remove last wash almost completely with the pipette.
- j. Prepare **Elution Buffer (EB)** and set the thermomixer to 65°C.
- k. Add 100 µl of **Elution Buffer** to each sample. Incubate for 10 minutes at 65°C in the thermomixer.
- l. Put the supernatant in a new tube and repeat the step k to reach a 200 µl final volume.
- m. Take 50 µl of the exceeded chromatin from the pre-clearing as a 5% input. Add 150 µl **Elution Buffer** to reach a 200 µl final volume.
- n. Incubate samples and inputs at 65°C overnight to de-crosslink.

Day 4

- a. Add 1 µl of Proteinase K to reach a 50 µg/ml final concentration. Incubate at 42°C for 1 hour.
- b. Elute samples twice with 30 µl of TE/EB/water until a final volume of 60 µl is reached.

CHROMATIN IMMUNOPRECIPITATION

www.ptglab.com

Buffer Solutions

1X PBS

10 mM Na₂HPO₄

1.8 mM KH₂PO₄

137 mM NaCl

2.7 mM KCl

Adjust pH to 7.4

Add ddH₂O to the final volume

Cell Lysis Buffer

5 mM HEPES

85 mM KCl

0.5% NP40, pH 8.0

Nuclear Lysis Buffer

50 mM Tris•HCl

10 mM EDTA

1% SDS, pH 8.1

Nuclear Lysis Buffer

0.1% SDS (protein interaction depend)

1.1% Triton X-100

1.2 mM EDTA

165 mM NaCl

16.7 mM Tris•HCl, pH 8.1

Low Salt Buffer

Tris•HCl 50 mM, pH 8.0

150 mM NaCl

0.1% SDS

1% NP40

1 mM EDTA

0.5% Deoxycholate Na

CHROMATIN IMMUNOPRECIPITATION

www.ptglab.com

Buffer Solutions

High Salt Buffer

Tris•HCl 50 mM, pH 8.0

500 mM NaCl

0.1% SDS

1% NP40

1 mM EDTA

0.5% Deoxycholate Na

LiCl Buffer

Tris•HCl 50 mM, pH 8.0

250 mM LiCl

0.1% SDS

1% NP40

1 mM EDTA

0.5% Deoxycholate Na

TE Buffer

Tris•HCl 10 mM, pH 8.0

0.25 mM EDTA

Elution Buffer

100 mM NaHCO₃

1% SDS