

IMMUNOHISTOCHEMISTRY USING SLIDE-MOUNTED PARAFFIN SECTIONS

www.ptglab.com

All steps in the following protocol are carried out at room temperature unless stated otherwise. Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

-
- 1. Deparaffinizing and rehydration:**
 - Immerse slides in xylene for 10 minutes. Repeat this step again in fresh xylene for 10 minutes. (If required, repeat a third time in fresh xylene for another 10 minutes.)
 - Rehydrate sections by sequentially incubating with 100%, 95%, 80% and 60% ethanol for 5 minutes each.
 - Rinse sections with distilled water three times for 3 minutes each.

 - 2. Antigen retrieval (optional):**
 - Transfer slides to a microwave-proof container and cover with **citrate buffer** or **Tris-EDTA (TE) buffer**.
 - Heat in the microwave on medium power for 10 minutes.
 - Allow slides to cool in the **citrate buffer** or **Tris-EDTA (TE) buffer** for approximately 35 minutes.

 - 3. Proteintech antibody incubation:**
 - Rinse slides three times with 1x **TBST** for 3 minutes each.
 - Incubate slides with 3% H₂O₂ solution (diluted in distilled water) for 10 minutes to quench endogenous peroxidase activity.
 - Rinse slides three times with 1x **TBST** for 3 minutes each, then rinse slides three times with distilled water for 3 minutes each.
 - Prepare 5% normal blocking serum in 1x **TBST**. The serum should be derived from the same species in which the secondary antibody was raised. Block the sections for 1 hour. (Alternatively, use 5% BSA in 1x **TBST** for blocking if the corresponding serum is not available.)
 - Incubate sections with primary antibody diluted in 1x **TBST** for 1 hour, or overnight at 4°C; the optimal antibody dilution ratio should be pre-determined by experimentation. Set up negative controls by omitting the primary antibody incubation step for one slide per each experimental condition.
 - Following primary antibody incubation rinse slides three times with 1x **TBST** for 3 minutes each.

IMMUNOHISTOCHEMISTRY USING SLIDE-MOUNTED PARAFFIN SECTIONS

www.ptglab.com

All steps in the following protocol are carried out at room temperature unless stated otherwise. Recipes for all solutions (highlighted) in **bold** are included at the end of the protocol.

-
- 4.** **Signal detection:**
Proteintech routinely uses EnVision Kit reagents (Dako) for this step.
- Apply sufficient peroxidase labeled polymer and incubate for 30 minutes.
 - Rinse slides three times with 1x **TBST** for 3 minutes each.
 - Prepare an appropriate volume of substrate solution prior to use by mixing one drop of Liquid DAB plus chromogen immediately with 1 ml of substrate buffer. Apply the substrate carefully and incubate for 5–10 minutes till a brown color develops.
 - Rinse sections gently with sufficient distilled water.
-
- 5.** **Hematoxylin counterstaining (optional):**
- To stain nuclei, immerse slides in a bath of hematoxylin for 3 minutes.
 - Rinse slides gently with distilled water.
 - Transfer slides into a 1% HCl, 99% ethanol solution for 10 seconds; transfer to distilled water immediately.
-
- 6.** **Dehydration and mounting:**
- Immerse slides sequentially into 60%, 80%, 95% and 100% ethanol baths for 5 minutes each.
 - Immerse slides in xylene for 5 minutes. Repeat this step again in fresh xylene for 5 minutes.
 - Mount the section with sufficient mounting media and cover with a cover slip. Air-dry in a well-ventilated area (e.g. fume hood).

IMMUNOHISTOCHEMISTRY USING SLIDE-MOUNTED PARAFFIN SECTIONS

www.ptglab.com

Solutions

Citrate buffer	For 1000 ml	1x TBS	For 1000 ml
10 mM Trisodium citrate+2H ₂ O	2.9 g	20 mM Tris-base	2.4 g
1.9 mM Citric acid+H ₂ O	0.4 g	150 mM NaCl	8.7 g
Adjust pH to 6.0		Adjust pH to 7.6	
Add ddH ₂ O to 1000 ml		Add ddH ₂ O to 1000 ml	

Tris-EDTA (TE) buffer	For 1000 ml	1x TBST	For 1000 ml
10 mM Tris-base	1.21 g	1x TBS	999 ml
1 mM EDTA C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ +2H ₂ O	0.372 g	Tween-20	1 ml
Adjust pH to 9.0			
Add ddH ₂ O to 1000 ml			