

Human CD19 Magnetic Beads Kit

Kit Contents (Cat #:KMS004)

- 1x biotinylated human
 CD19 antibody (Cat #:MS65197)
- 1x streptavidin magnetic beads (Cat #:MS001)



Protocol

- 1. Take the cells of interest, wash and re-suspend in cell separation buffer PBS, 0.1% BSA, 2mM EDTA, pH 7.4 (100μL for every 10⁷ cells).
- 2. Incubate cells with $10\mu L$ biotin conjugated antibody at 4°C for 30 minutes.
- 3. Wash cells with 2mL PBS and re-suspend in $100\mu L$ cell separation buffer.
- 4. Add $10\mu L$ of streptavidin magnetic beads in each tube and incubate at 4°C for 30 minutes.
- 5. After incubation add 2mL PBS in the tube and put the tube on magnet rack for 10 minutes.
- 6. Gently remove supernatant, avoiding contact with the cells bound to magnetic beads.
- 7. The supernatant contains the depleted cells, the enriched cells remain in the tube.
- 8. Remove tube from magnet, re-suspend cells in 2mL PBS and wash.
- 9. Now your cells are ready for further analysis.
- 10. If required, repeat steps 2-8 on the enriched cells for better results.